

DAPI Staining of Bacteria

Updated March 2017 – Rachel Parsons

This experiment is to enumerate bacteria using a microscope. To do this you need to get the bacteria onto a flat surface. A filter provides this surface. We use a 0.2 μ M filter that allows the water and viruses to pass through the membrane but retains the bacteria and larger plankton onto the filter so they can be counted.

Equipment

Filtration rig – usually 4-5 place
Filter bases
Filter Towers
Vacuum Flask
Pump
Backing Filters – Supor 0.8 μ M 25mm
Filters – PC 0.2 μ M 25mm
Pipet – 5-10ml
Pipet Tips 5-10ml
Pipet 1ml
Pipet tips 1ml
Forceps
Slides
Cover slips
High Viscosity Oil
Ultra fine sharpie
50ml falcon tubes
0.2 μ M filtered formalin
Glass waste container
Fixative Waste container – Haz Mat labeled.
Epifluorescent Microscope with UV light filter cubes.


Reagents

DAPI 5 μ g/ml
Stock is 1ml of 200 μ g/ml. Dilute into 40mls 0.2 μ M filtered QW using 50ml falcon tube.

Irgalan Black for staining the filters.
200mg of Irgalan Black dye into 2% Acetic Acid (2mls Acetic Acid into 98mls of 0.2 μ M filtered QW)

Q Water to rinse the filters – 0.2 μ M filtered

Filtration set-up:



Start of filtration:

- Label slides in the following way;

Sample ID (Cruise)	→	B183
Depth(m)		0m 10m
		A B
mls filtered/tower diameter	→	10mls 10mls

- Take samples out of the fridge and shake them well. Allow them to come to room temperature. The samples have been killed with formalin (2%) to stop bacterial replication.
- Mount 0.8um backing filter (GN4 Metrice) followed by blackened 0.2um PC filter on all 4 filter towers. Leave all valves open.
- Pipet 10 mls (BATS surface) or 3-5mL (Bermuda Inshore) of sample into the filter tower using the DEAD pipets in the fume hood.
- Start pump (vacuum on 5-8 with all valves closed).
- Let filter until all sample on all towers has gone through.
- While samples are filtering prepare slides for mounting by putting a drop of immersion oil on the slide and spreading it across the slide with a cover slip.
- Stop filtration by stopping the pump using the foot control.
- Add approximately 0.5-1ml DAPI solution to each sample using the DEAD pipet found in the fume hood. DAPI is a stain that binds to DNA and fluoresces under UV light so the cells can be seen under the microscope.
- Let stain for 3 minutes (use the timer provided).
- Turn on pump to draw down DAPI.
- As soon as DAPI goes through, take the filter off the tower while the pump and vacuum are still on. Mount two filters per slide, add a drop of immersion oil on the filter, and cover with a cover slip.
- Store prepared filters in slide-box at room-temp. until finished preparing the samples.
- Put slide box in DEAD fridge until you examine them under the microscope.